

Rocker Chemically Resistant Vacuum Pump



All wetted parts of Chemker series pump are made of PTFE material which can resist most of corrosive gases.

Model	Rocker 300C	Chemker 300	Chemker 400	Chemker 410
Max Vacuum (mmHg)	-630	-670	-670	-750
Flow Rate (l/min)	20	18	33	18
Dimensions (LxWxH) cm	28.5 x 13.5 x 20.4	20.3x11.2x20.9	26.7 x 11.2 x 18.1	26.2 x 14.8 x18.2
Weight (Kg)	4.1	5.0	7.0	7.0
Noise Level (dB)	50	50	50	50

Serving Science Since 1975

www.severnsaleslabequip.com



+44 (0) 117 9608858
sales@severnsaleslabequip.com