

Rocker Oil Free Vacuum Pumps



Rocker series pumps are driven by piston, without the need of lubricant, regular oil changes and maintenance; with no oil pollution.

Model	Rocker 300	Rocker 400	Rocker 410	Rocker 800
Max Vacuum (mmHg)	-680	-680	-735	-670
Flow Rate (l/min)	20	34	20	80
Dimensions (LxWxH) cm	26.8 x 13.5 x 20.4	31.0 x 13.5 x 20.4	31.0 x 13.5 x 20.4	34.5 x 17.5 x 20.0
Weight (Kg)	4.1	5.2	5.2	7.6
Noise Level (dB)	50	5.2	52	68

Serving Science Since 1975

www.severnsaleslabequip.com



+44 (0) 117 9608858
sales@severnsaleslabequip.com