

47 mm Stainless Steel Filter Holder



Features

- Spin-lock connection facilitates fast and stable installation without clamp.
- Stainless steel funnel offers excellent chemical resistance and can be fastly sterilized by flame.
- Funnel lid is optional.

Specifications

Material	100ml	300ml	500ml
Funnel	SS316	SS316	SS316
Support base	SS316	SS316	SS316
Membrane support	SS316	SS316	SS316
Funnel lid	SS316	SS316	SS316

Specifications

Funnel capacity	100 ml	300 ml	500 ml
Filter diameter	47 mm / 50 mm	47 mm / 50 mm	47 mm / 50 mm
Effective filtration area	9.6 cm ²	9.6 cm ²	9.6 cm ²
SS membrane pore size	70 um	70 um	70 um

Ordering Information

100 ml, 47 mm Stainless Steel Filter Holder 167103-10

1	Stainless steel funnel, 100 ml	167103-61
2	Stainless steel support screen	167100-31
3	Stainless steel support base	167103-60
4	Silicone stopper (No.8)	167110-16

300 ml, 47 mm Stainless Steel Filter Holder 167103-30

1	Stainless steel funnel, 300 ml	167103-63
2	Stainless steel support screen	167100-31
3	Stainless steel support base	167103-60
4	Silicone stopper (No.8)	167110-16

500 ml, 47 mm Stainless Steel Filter Holder 167103-50

1	Stainless steel funnel, 500 ml	167103-65
2	Stainless steel support screen	167100-31
3	Stainless steel support base	167103-60
4	Silicone stopper (No.8)	167110-16

Stainless Steel Lid

100 ml Stainless steel lid	167103-11
300 ml Stainless steel lid	167103-13
500 ml Stainless steel lid	167103-15

