



ULTF130 -80°C Ultra Low Temperature Sparkfree Chest Freezer

Compact ultra low chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and castors for easy manoeuvrability.



Featuring an energy saving Danfoss compressor this freezer offers little temperature variation across the cabinet, making it ideal for the storage of samples which demand tight temperature tolerances.

Bespoke racking and cryoboxes can be provided as an optional extra and to prolong the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Part Number	ULTF130
Capacity (Litres)	130
Style	Chest Freezer
Exterior Dimensions (HxWxD)	920 x 720 x 735mm
Interior Dimensions (HxWxD)	650 x 450 x 520mm
Weight	56kg
Operating Temperature	-65°C to -85°C
Interior Construction	PVC Coated Steel
Exterior Construction	Painted Steel
Defrost	Manual
Access Port	✓
Sparkfree Interior	✓
Open Door Alarm	x
High Temperature Alarm	✓
Low Temperature Alarm	✓
Mains Failure Alarm	✓
Remote Alarm Contacts	✓
Digital Temperature Display	✓
Key Operated Lid Lock	✓
Internal Polystyrene Lids	✓
Castors	4
Vacuum Release Port	✓
Noise Level	≤51dba
Energy Consumption	5.8 kWh/24h
Eco Friendly Refrigerant	✓